## BLEACH MASK VD

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## BLEACH MASK VD

### A-Silicone for Gingiva Protection during Bleaching

BLEACH MASK VD is a fast setting addition curing vinyl silicone for isolating the gingiva adjacent to teeth being whitened. It shows excellent flow properties for precise positioning and for sealing interproximal spaces.

Delivered in MINIMIX-syringes, BLEACH MASK VD can be applied directly and bubble free onto the gingiva. After tooth whitening it is easy to remove in one or a few large pieces.

### Preparing of the MINIMIX-syringe

Take off the turn cap of the syringe and replace it by a mixing cannula (do not use it again!).

#### Note:

**Do not re-use the syringe cap.** Closing the syringe by turning on the cap between use could cause cross contamination of components and thus premature setting.

### Application of BLEACH MASK VD

Ensure that teeth to be covered are cleaned thoroughly. Apply BLEACH MASK VD to build a 4-6 mm wide strip on the gingiva and seal the interproximal spaces. Overlap BLEACH MASK VD approximately 0.5 mm onto dry enamel.

Extend BLEACH MASK VD one tooth beyond the last tooth to be bleached. Visually check that all gingival tissues at BLEACH MASK VD margin are covered and seal is established.

BLEACH MASK VD has a working time of ca. 1:30 minutes (working time at 23°C (74°F)). The setting time in mouth is 2:00 minutes. After setting has finished, the bleaching gel is applied according to user instructions.

After bleaching procedure is finished remove the BLEACH MASK VD and polish the teeth.

### Storage

Do not store above 78°F (25°C). Do not use after expiration date.

#### Additional notes:

Do not use Latex gloves, they may interfere with the setting reaction of the impression material.