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## Instructions for use

Bleach 'n Smile bleaching gel is intended solely for bleaching in the dental practice.

### Indications:

- whitening devitalised teeth
- whitening teeth discoloured by disease, as an alternative to a crown or veneer
- treatment of teeth discoloured by medication (e.g. tetracycline)

### Contraindications:

Exposed dentine, defective filling margins, hypersensitive patients. Some patients experience a slight sensitisation of the teeth during or shortly after treatment, but this normally lasts only a short time. **Treatment should be terminated immediately if the patient experiences excessive pain during bleaching.**

### Preparation:

Determine the shade using the Vita® shade guide. (Foto 1):  
For your assessment, rearrange the Vita shade-guide as follows:  
Values arranged according to the shade: A1 A2 A3 A3,5 A4 B1 B2 B3 B4 C1 C2 C3 C4 D2 D3 D4  
Arranged according to brightness: B1 A1 B2 D2 A2 C1 C2 D3 A3 D4 B3 A3,5 B4 C3 A4 C4  
B2 to B1 = increase by 2, A3 to A2= increase by 4  
Clean the patient's teeth with an Air-Flow or prophylaxis paste. Insert a cheek or lip retractor or similar instrument. Ensure that the gingiva is protected., e.g. by using Bleach 'n Smile Dental Dam or placing a rubber dam. Do not anaesthetise the patient. Bleach 'n Smile bleaching gel and activator are caustic. Wear gloves and protective glasses when using them. It is essential to avoid contact with the eyes, skin and oral tissue.

### Mixing the gel:

Manually twist off the tip of the plastic container (content: red activator fluid). Add the liquid in drops to the pot with the bleaching gel. (Foto 2) Then mix the gel thoroughly for approx. 1 min. with the spatula provided, until it is a homogenous red colour. (Foto 3) The gel is now ready to use and can be applied to the teeth to be bleached using the spatula and/or the brush supplied. Do not mix partial quantities, but always the complete package. Mixing partial quantities could cause discolorations or insufficient bleaching results.

### Application:

The gel should be applied approx. 1-2 mm thick to the labial surfaces. (Foto 4) A standard plasma light-curing unit can be used to increase the effect of the gel. The patient and operator should wear protective glasses during light curing. Follow the manufacturer's instructions when light-curing. Overheating must be avoided.

Allow the gel to take effect for approx. 10 min., then remove it with suction and wipe off any residue. If a greater degree of bleaching is required, the procedure can be repeated twice more until the desired effect is achieved. Ensure beforehand that there is still sufficient dental dam. (Foto 5-6)

Then check the shade again. (Foto 7)

Note: any existing composite fillings as well as porcelain or metal restorations are not bleached. Fillings should therefore be renewed 2 weeks after treatment.

### Devitalized teeth:

The gel is also applied via a lingual access cavity. Fill the gel into an expendable syringe and apply it through a tip. To do so, remove the root canal filling to within approx. 2 mm of the cementoenamel junction. Then cover the gel with glass ionomer cement, ensuring that the margin is sealed. Bleach the tooth as described previously.

### Finishing:

After the final bleaching procedure, first remove the gel with suction, wipe it away and then thoroughly rinse off with water. Then dry the surface of the teeth and polish them with polishing paste. Minerals should then be replaced using a fluoride solution or gel. These can also be used for relieving pain if sensitisation occurs. After treatment, patients should not eat, drink or smoke for 1 hour and in the next few days they should avoid strongly coloured substances such as red wine.

### Shelf life and storage:

Bleach 'n Smile gel should not be used after the expiry date, as an adequate bleaching effect can no longer be guaranteed. Any mixed Bleach 'n Smile gel (bleaching gel and activator) should be used on the same day and certainly no later than the following day. Store Bleach 'n Smile gel in the refrigerator at 3°C- 9°C.

**Product contains 35% hydrogen peroxide, potassium nitrate, sodium fluoride and glycol. It should not be used in cases known to be allergic or intolerant to either of these substances.**

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Schütz Dental GmbH